

## ***DNAnimal Kits***

Animal Species Testing Made Easy





## DNAnimal Kits for Your Lab

The DNAnimal kit line enables you to identify and differentiate animal species in your sample DNA.

DNAnimal kits are available for most domestic animal species as well as for some specific applications like halal testing or animal species groups i.e. ruminants for feed testing or mammals and birds i.e. for testing of vegan food.

Eurofins GeneScan is the central production unit for kits and reagents for the labs of the Eurofins group as well as customers around the world – benefit from kits derived by the experience and quality of a market leader in molecular biological testing!

## DNA Extraction Kits

For DNA extraction, we offer a number of different kits for different matrices. Choose the optimal method for your sample matrix.

## Versatile Technology

All kits are validated on ABI7500, Bio-Rad CFX96, Agilent MX3005P and AriaMX. Some are also validated on Roche LC480 and ABI7900. Most kits contain an IPC (internal positive control) which protects against false-negative results.

All kits detect the species in question with a sensitivity of 0.01%.

## Our Offer

Our portfolio is constantly growing. Currently, the following kits are available:

- DNAnimal Screen Halal (horse, donkey and pig)
- DNAnimal Screen Tetrapod (mammals and birds)
- DNAnimal Screen Ruminant
- DNAnimal Ident Beef
- DNAnimal Ident Pork
- DNAnimal Ident Horse
- DNAnimal Ident Water Buffalo
- DNAnimal Ident Sheep
- DNAnimal Ident Goat
- DNAnimal Ident Chicken
- DNAnimal Ident Duck
- DNAnimal Ident Goose
- DNAnimal Ident Turkey

Furthermore, we offer a broad range of animal species DNA as positive control DNA and meat testing standards with defined percentages of i.e. horse or pig.