

SureID[®] 21G Human STR Identification Kit

- **High compatibility**

Direct sample amplification for Whatman blood/saliva - FTA cards, normal blood/saliva cards, oral swabs, purified DNA samples, etc.

- **Optimized peak balance**

Optimized peak balance of STRs, especially heterozygous alleles.


- **Optimized reaction system**

Minimize amplification time to 90 min.

- **High stability**


Components can be stored for longer.





- Allelic Ladder Mix
- Control DNA 9947A
- Normal blood card
- Normal saliva card
- Filter paper
- Blood-FTA card
- Saliva-FTA card
- Gauze
- Hair

Suited for various samples



Filter paper sample

Components

Pre-amplification Reagents	PCR Master Mix
	SureID® 21G Primer Mix
	Control DNA 9947A
	DNase/RNase-Free H ₂ O
Post-amplification Reagents	SureID® 21G Allelic Ladder Mix
	SIZE-500

Order Information

Product	Size	Item #
SureID® 21G Human STR Identification Kit	25 µL × 200 reactions/box	1060104

Forensic products are NOT available in USA market.